

## Product datasheet for **TP726099**

### **POLR3A Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human RNA polymerase III subunit C1
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	K1071-P1370
<b>Tag:</b>	C-His
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	MADE TO ORDER   Recombinant Human RNA polymerase III subunit C1 is produced by E.coli expression system and the target gene encoding K1071-P1370 is expressed with a His tag at the C-terminus.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	11128
<b>UniProt ID:</b>	<u><a href="#">O14802</a></u>



[View online »](#)