

Product datasheet for **TP726101**

DARS2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Aspartate--tRNA ligase, mitochondrial
Species:	Human
Expression cDNA Clone or AA Sequence:	A406-L626
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Aspartate--tRNA ligase, mitochondrial is produced by E.coli expression system and the target gene encoding A406-L626 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	55157
UniProt ID:	Q6PI48



[View online »](#)