

Product datasheet for **TP726102**

SSH3BP1 (ABI1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Abl interactor 1
Species:	Human
Expression cDNA Clone or AA Sequence:	G380-D508
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Abl interactor 1 is produced by E.coli expression system and the target gene encoding G380-D508 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10006
UniProt ID:	Q8IZP0
Synonyms:	Abelson interactor 1 (Abi-1); Abl-binding protein 4 (AbIBP4); Eps8 SH3 domain-binding protein (Eps8-binding protein); Nap1-binding protein (Nap1BP); Spectrin SH3 domain-binding protein 1; e3B1



[View online »](#)