

## Product datasheet for **TP726111**

### RPL15 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human 60S ribosomal protein L15
Species:	Human
Expression cDNA Clone or AA Sequence:	G2-R204
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human 60S ribosomal protein L15 is produced by E.coli expression system and the target gene encoding G2-R204 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	6138
UniProt ID:	<a href="#">P61313</a>
Synonyms:	Large ribosomal subunit protein eL15
Protein Pathways:	Ribosome



[View online »](#)