

## Product datasheet for **TP726112**

### STK25 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Serine/threonine-protein kinase 25
Species:	Human
Expression cDNA Clone or AA Sequence:	N230-R426
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Serine/threonine-protein kinase 25 is produced by E.coli expression system and the target gene encoding N230-R426 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10494
UniProt ID:	<a href="#">O00506</a>



[View online »](#)