

Product datasheet for **TP726115**

ornithine aminotransferase (OAT) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ornithine aminotransferase, mitochondrial
Species:	Human
Expression cDNA Clone or AA Sequence:	A173-E430
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Ornithine aminotransferase, mitochondrial is produced by E.coli expression system and the target gene encoding A173-E430 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	4942
UniProt ID:	P04181
Synonyms:	Ornithine delta-aminotransferase; Ornithine--oxo-acid aminotransferase
Protein Families:	Druggable Genome
Protein Pathways:	Arginine and proline metabolism, Metabolic pathways



[View online »](#)