

Product datasheet for **TP726118**

FGF14 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human FGF-14
Species:	Human
Expression cDNA Clone or AA Sequence:	D65-T247
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human FGF-14 is produced by E.coli expression system and the target gene encoding D65-T247 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	2259
UniProt ID:	Q92915
Synonyms:	Fibroblast growth factor 14; Fibroblast growth factor homologous factor 4 (FHF-4)
Protein Families:	Secreted Protein
Protein Pathways:	MAPK signaling pathway, Melanoma, Pathways in cancer, Regulation of actin cytoskeleton



[View online »](#)