

## Product datasheet for **TP726123**

### WDR62 Human Recombinant Protein

#### Product data:

|                                       |  |
|---------------------------------------|--|
| Product Type:                         | Recombinant Proteins   |
| Description:                          | Recombinant Human WD repeat-containing protein 62  |
| Species:                              | Human  |
| Expression cDNA Clone or AA Sequence: | N750-D900  |
| Tag:                                  | C-His  |
| Buffer:                               | Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4  |
| Note:                                 | MADE TO ORDER   Recombinant Human WD repeat-containing protein 62 is produced by E.coli expression system and the target gene encoding N750-D900 is expressed with a His tag at the C-terminus.              |
| Storage:                              | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles. |
| Stability:                            | 12 months from date of despatch  |
| Locus ID:                             | 284403   |
| UniProt ID:                           | <a href="#">O43379</a>   |
| Synonyms:                             | WDR62 Protein  |



[View online »](#)