

Product datasheet for **TP726124**

CDX1 Human Recombinant Protein

Product data:

| | |
|---------------------------------------|--|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human Homeobox protein CDX-1 |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | R137-P265 |
| Tag: | C-His |
| Buffer: | Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | MADE TO ORDER Recombinant Human Homeobox protein CDX-1 is produced by E.coli expression system and the target gene encoding R137-P265 is expressed with a His tag at the C-terminus. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 1044 |
| UniProt ID: | P47902 |



[View online »](#)