

Product datasheet for TP726125

CYP2A13 Human Recombinant Protein

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Recombinant Proteins
Description:	Recombinant Human Cytochrome P450 2A13/CYP2A13 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	P350-R494
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Cytochrome P450 2A13/CYP2A13 Protein is produced by E.coli expression system and the target gene encoding P350-R494 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	1553
UniProt ID:	<u>Q16696</u>



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US