

Product datasheet for **TP726128**

ALDH6A1 Human Recombinant Protein

Product data:

| | |
|--|---|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human MMSDH |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | V326-T495 |
| Tag: | C-His |
| Buffer: | Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | MADE TO ORDER Recombinant Human MMSDH is produced by E.coli expression system and the target gene encoding V326-T495 is expressed with a His tag at the C-terminus. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 4329 |
| UniProt ID: | <u>Q02252</u> |



[View online »](#)