

## Product datasheet for **TP726133**

### PPP2R2B Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human PPP2R2B Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	V181-E380
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human PPP2R2B Protein is produced by E.coli expression system and the target gene encoding V181-E380 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	5521
UniProt ID:	<a href="#">Q00005</a>
Synonyms:	Serine/threonine-protein phosphatase 2A 55 kDa regulatory subunit B beta isoform;PP2A subunit B isoform B55-beta; PP2A subunit B isoform PR55-beta; PP2A subunit B isoform R2-beta; PP2A subunit B isoform beta
Protein Families:	Druggable Genome, Phosphatase
Protein Pathways:	Tight junction



[View online »](#)