

Product datasheet for **TP726137**

BICRA Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human BICRA Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	H1331-E1499
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human BICRA Protein is produced by E.coli expression system and the target gene encoding H1331-E1499 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	29998
UniProt ID:	Q9NZM4
Synonyms:	BRD4-interacting chromatin-remodeling complex-associated protein; Glioma tumor suppressor candidate region gene 1 protein
Protein Families:	Druggable Genome



[View online »](#)