

Product datasheet for **TP726145**

Acetoacetyl CoA synthetase (AACS) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Acetoacetyl-CoA synthetase
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-I239
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Acetoacetyl-CoA synthetase is produced by E.coli expression system and the target gene encoding S2-I239 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	65985
UniProt ID:	Q86V21



[View online »](#)