

Product datasheet for **TP726147**

CRPPA Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human CRPPA Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	R22-A323
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human CRPPA Protein is produced by E.coli expression system and the target gene encoding R22-A323 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	729920
UniProt ID:	A4D126
Synonyms:	D-ribitol-5-phosphate cytidyltransferase; 2-C-methyl-D-erythritol 4-phosphate cytidyltransferase-like protein; Isoprenoid synthase domain-containing protein



[View online »](#)