

Product datasheet for **TP726150**

ARSL Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Arylsulfatase L
Species:	Human
Expression cDNA Clone or AA Sequence:	E308-L453
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Arylsulfatase L is produced by E.coli expression system and the target gene encoding E308-L453 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	415
UniProt ID:	P51690
Synonyms:	Arylsulfatase E (ASE)
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)