

Product datasheet for **TP726151**

WDPCP Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human FRTZ
Species:	Human
Expression cDNA Clone or AA Sequence:	R533-K739
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human FRTZ is produced by E.coli expression system and the target gene encoding R533-K739 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	51057
UniProt ID:	O95876
Synonyms:	WD repeat-containing and planar cell polarity effector protein fritz homolog; Bardet-Biedl syndrome 15 protein; WD repeat-containing and planar cell polarity effector protein



[View online »](#)