

## Product datasheet for **TP726155**

### ADH6 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Alcohol dehydrogenase 6
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-D264
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Alcohol dehydrogenase 6 is produced by E.coli expression system and the target gene encoding S2-D264 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	130
UniProt ID:	<a href="#">P28332</a>



[View online »](#)