

Product datasheet for **TP726159**

CCT3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human TCP-1-gamma
Species:	Human
Expression cDNA Clone or AA Sequence:	E246-E545
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human TCP-1-gamma is produced by E.coli expression system and the target gene encoding E246-E545 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	7203
UniProt ID:	P49368
Synonyms:	T-complex protein 1 subunit gamma; CCT-gamma; hTRiC5



[View online »](#)