

## Product datasheet for **TP726161**

### TCP1 theta (CCT8) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human TCP-1-theta
Species:	Human
Expression cDNA Clone or AA Sequence:	F197-D548
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human TCP-1-theta is produced by E.coli expression system and the target gene encoding F197-D548 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10694
UniProt ID:	<a href="#">P50990</a>
Synonyms:	T-complex protein 1 subunit theta; CCT-theta; Chaperonin containing T-complex polypeptide 1 subunit 8; Renal carcinoma antigen NY-REN-15
Protein Families:	Druggable Genome



[View online »](#)