

Product datasheet for **TP726165**

LIMS1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human LIMS1
Species:	Human
Expression cDNA Clone or AA Sequence:	I83-V204
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human LIM and senescent cell antigen-like-containing domain protein 1 is produced by E.coli expression system and the target gene encoding I83-V204 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	3987
UniProt ID:	<u>P48059</u>



[View online »](#)