

Product datasheet for **TP726166**

DPP8 Human Recombinant Protein

Product data:

| | |
|---------------------------------------|--|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human DP8 |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | A18-L135 |
| Tag: | C-His |
| Buffer: | Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | MADE TO ORDER Recombinant Human DP8 is produced by E.coli expression system and the target gene encoding A18-L135 is expressed with a His tag at the C-terminus. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 54878 |
| UniProt ID: | Q6V1X1 |



[View online »](#)