

Product datasheet for **TP726167**

CPNE1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Copine-1
Species:	Human
Expression cDNA Clone or AA Sequence:	K25-S220
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Copine-1 is produced by E.coli expression system and the target gene encoding K25-S220 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	8904
UniProt ID:	Q99829



[View online »](#)