

Product datasheet for **TP726173**

Acetyl CoA synthetase (ACSS2) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Acetyl-coenzyme A synthetase, cytoplasmic/ ACSS2 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	P521-S690
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Acetyl-coenzyme A synthetase, cytoplasmic/ ACSS2 Protein is produced by E.coli expression system and the target gene encoding P521-S690 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	55902
UniProt ID:	Q9NR19



[View online »](#)