

Product datasheet for **TP726177**

Triosephosphate isomerase (TPI1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human TIM
Species:	Human
Expression cDNA Clone or AA Sequence:	G129-F241
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human TIM is produced by E.coli expression system and the target gene encoding G129-F241 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	7167
UniProt ID:	P60174



[View online »](#)