

Product datasheet for **TP726180**

SIVA (SIVA1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Apoptosis regulatory protein Siva
Species:	Human
Expression cDNA Clone or AA Sequence:	P2-A110
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Apoptosis regulatory protein Siva is produced by E.coli expression system and the target gene encoding P2-A110 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10572
UniProt ID:	O15304
Synonyms:	CD27-binding protein (CD27BP)
Protein Families:	Druggable Genome



[View online »](#)