

Product datasheet for **TP726181**

MBNL2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Muscleblind-like protein 2
Species:	Human
Expression cDNA Clone or AA Sequence:	P160-A260
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Muscleblind-like protein 2 is produced by E.coli expression system and the target gene encoding P160-A260 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	10150
UniProt ID:	Q5VZF2



[View online »](#)