

## Product datasheet for **TP726182**

### MICAL1 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human [F-actin]-monooxygenase MICAL1
Species:	Human
Expression cDNA Clone or AA Sequence:	H718-G1067
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human [F-actin]-monooxygenase MICAL1 is produced by E.coli expression system and the target gene encoding H718-G1067 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	64780
UniProt ID:	<a href="#">Q8TDZ2</a>



[View online »](#)