

## Product datasheet for **TP726184**

### USP45 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Ubiquitin carboxyl-terminal hydrolase 45
Species:	Human
Expression cDNA Clone or AA Sequence:	N323-P536
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Ubiquitin carboxyl-terminal hydrolase 45 is produced by E.coli expression system and the target gene encoding N323-P536 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	85015
UniProt ID:	<a href="#">Q70EL2</a>
Synonyms:	Deubiquitinating enzyme 45; Ubiquitin thioesterase 45; Ubiquitin-specific-processing protease 45
Protein Families:	Druggable Genome, Protease



[View online »](#)