

Product datasheet for **TP726186**

CRTC3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human CREB-regulated transcription coactivator 3
Species:	Human
Expression cDNA Clone or AA Sequence:	S425-N574
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human CREB-regulated transcription coactivator 3 is produced by E.coli expression system and the target gene encoding S425-N574 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	64784
UniProt ID:	Q6UUV7
Synonyms:	Transducer of regulated cAMP response element-binding protein 3 (TORC-3; Transducer of CREB protein 3)



[View online »](#)