

Product datasheet for **TP726188**

ZNF313 (RNF114) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human E3 ubiquitin-protein ligase RNF114
Species:	Human
Expression cDNA Clone or AA Sequence:	E119-Q228
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human E3 ubiquitin-protein ligase RNF114 is produced by E.coli expression system and the target gene encoding E119-Q228 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	55905
UniProt ID:	Q9Y508



[View online »](#)