

Product datasheet for **TP726188**

ZNF313 (RNF114) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human E3 ubiquitin-protein ligase RNF114
Species:	Human
Expression cDNA Clone or AA Sequence:	E119-Q228
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human E3 ubiquitin-protein ligase RNF114 is produced by E.coli expression system and the target gene encoding E119-Q228 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	55905
UniProt ID:	<u>Q9Y508</u>

[View online »](#)