

## Product datasheet for **TP726189**

### COPS4 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human SGN4
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-L230
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human SGN4 is produced by E.coli expression system and the target gene encoding A2-L230 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	51138
UniProt ID:	<a href="#">Q9BT78</a>
Synonyms:	COP9 signalosome complex subunit 4; Signalosome subunit 4; JAB1-containing signalosome subunit 4



[View online »](#)