

Product datasheet for **TP726191**

MAB21L2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Protein mab-21-like 2
Species:	Human
Expression cDNA Clone or AA Sequence:	V200-L359
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Protein mab-21-like 2 is produced by E.coli expression system and the target gene encoding V200-L359 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10586
UniProt ID:	Q9Y586
Synonyms:	MAB21L2 Protein



[View online »](#)