

Product datasheet for **TP726198**

PWWP3A Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human PWWP3A
Species:	Human
Expression cDNA Clone or AA Sequence:	E100-L280
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human PWWP3A is produced by E.coli expression system and the target gene encoding E100-L280 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	84939
UniProt ID:	Q2TAK8
Synonyms:	PWWP domain-containing DNA repair factor 3A; Mutated melanoma-associated antigen 1 (MUM-1); PWWP domain-containing protein MUM1; Protein expandere



[View online »](#)