

## Product datasheet for **TP726199**

### **NAB2 Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human NGFI-A-binding protein 2
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	I130-L320
<b>Tag:</b>	C-His
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	MADE TO ORDER   Recombinant Human NGFI-A-binding protein 2 is produced by E.coli expression system and the target gene encoding I130-L320 is expressed with a His tag at the C-terminus.
<b>Storage:</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	4665
<b>UniProt ID:</b>	<a href="#">Q15742</a>
<b>Synonyms:</b>	EGR-1-binding protein 2; Melanoma-associated delayed early response protein (Protein MADER)
<b>Protein Families:</b>	Transcription Factors



[View online »](#)