

Product datasheet for **TP726201**

KRT6A Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Keratin, type II cytoskeletal 6A
Species:	Human
Expression cDNA Clone or AA Sequence:	E162-L176
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 5%Glycerol, 1mM DTT, pH7.4
Note:	MADE TO ORDER Recombinant Human Keratin, type II cytoskeletal 6A is produced by E.coli expression system and the target gene encoding E162-L176 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	3853
UniProt ID:	P02538



[View online »](#)