

Product datasheet for **TP726202**

MYBPC2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Fast MyBP-C
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-Q438
Buffer:	Supplied as a 0.2 um filtered solution of PBS, pH7.4
Note:	MADE TO ORDER Recombinant Human Fast MyBP-C is produced by E.coli expression system and the target gene encoding M1-Q438 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	4606
UniProt ID:	Q14324
Synonyms:	Myosin-binding protein C, fast-type; C-protein, skeletal muscle fast isoform



[View online »](#)