

Product datasheet for **TP726204**

Osteocalcin (BGLAP) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Osteocalcin
Species:	Human
Expression cDNA Clone or AA Sequence:	K24-100V
Buffer:	Supplied as a 0.2 um filtered solution of PBS, pH7.4
Note:	MADE TO ORDER Recombinant Human Osteocalcin is produced by E.coli expression system and the target gene encoding K24-100V is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	632
UniProt ID:	P02818
Synonyms:	Bone Gla protein (BGP); Gamma-carboxyglutamic acid-containing protein
Protein Families:	ES Cell Differentiation/IPS, Secreted Protein



[View online »](#)