

Product datasheet for **TP726205**

myosin light chain 1 (MYL1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human MLC1F/MLC3F
Species:	Human
Expression cDNA Clone or AA Sequence:	E54-I194
Buffer:	Supplied as a 0.2 um filtered solution of PBS, pH7.4
Note:	MADE TO ORDER Recombinant Human MLC1F/MLC3F is produced by E.coli expression system and the target gene encoding E54-I194 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	4632
UniProt ID:	P05976
Synonyms:	Myosin light chain 1/3, skeletal muscle isoform; MLC1/MLC3; Myosin light chain alkali 1/2 (Myosin light chain A1/A2)



[View online »](#)