

Product datasheet for **TP726206**

Myosin Light Chain 2 (MYL2) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human MLC-2
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-D166
Buffer:	Supplied as a 0.2 um filtered solution of PBS, pH7.4
Note:	MADE TO ORDER Recombinant Human MLC-2 is produced by E.coli expression system and the target gene encoding M1-D166 is expressed with a His tag at the C-terminus.
Stability:	12 months from date of despatch
Locus ID:	4633
UniProt ID:	P10916



[View online »](#)