

Product datasheet for **TP726209**

VPS13B Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Intermembrane lipid transfer protein VPS13B
Species:	Human
Expression cDNA Clone or AA Sequence:	E64-Y412
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Intermembrane lipid transfer protein VPS13Bis produced by E. coli. The target gene encodingE64-Y412 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	157680
UniProt ID:	Q7Z7G8
Synonyms:	Intermembrane lipid transfer protein VPS13B; VPS13B; Cohen syndrome protein 1; Vacuolar protein sorting-associated protein 13B



[View online »](#)