

Product datasheet for **TP726210**

TUBA8 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Tubulin alpha-8 chain/TUBA8
Species:	Human
Expression cDNA Clone or AA Sequence:	Q31-F169
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Tubulin alpha-8 chain/TUBA8is produced by E. coli. The target gene encodingQ31-F169 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	51807
UniProt ID:	Q9NY65
Synonyms:	Tubulin alpha-8 chain; TUBA8; Alpha-tubulin 8; Tubulin alpha chain-like 2
Protein Families:	Druggable Genome
Protein Pathways:	Gap junction, Pathogenic Escherichia coli infection



[View online »](#)