

Product datasheet for **TP726213**

RABL4 (IFT27) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Intraflagellar transport protein 27 homolog/IFT27
Species:	Human
Expression cDNA Clone or AA Sequence:	T92-A186
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Intraflagellar transport protein 27 homolog/IFT27 is produced by E. coli. The target gene encoding T92-A186 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	11020
UniProt ID:	Q9BW83
Synonyms:	Intraflagellar transport protein 27 homolog; IFT27; Putative GTP-binding protein RAY-like; Rab-like protein 4
Protein Families:	Druggable Genome



[View online »](#)