

## Product datasheet for **TP726215**

### TRMT1 (TRMU) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Mitochondrial tRNA-specific 2-thiouridylase 1/TRMU
Species:	Human
Expression cDNA Clone or AA Sequence:	L22-W259
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Mitochondrial tRNA-specific 2-thiouridylase 1/TRMU is produced by E. coli. The target gene encoding L22-W259 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	55687
UniProt ID:	<a href="#">O75648</a>
Synonyms:	Mitochondrial tRNA-specific 2-thiouridylase 1; TRMU; MTO2 homolog



[View online »](#)