

Product datasheet for **TP726220**

CPNE3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Copine-3/CPNE3
Species:	Human
Expression cDNA Clone or AA Sequence:	A2-V266
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Copine-3/CPNE3is produced by E. coli. The target gene encodingA2-V266 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	8895
UniProt ID:	O75131
Synonyms:	Copine-3; CPNE3; Copine III



[View online »](#)