

Product datasheet for **TP726223**

PMP2 Human Recombinant Protein

Product data:

| | |
|---------------------------------------|--|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human Myelin P2 protein/PMP2 |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | S2-V132 |
| Tag: | C-His |
| Buffer: | Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | Recombinant Human Myelin P2 protein/PMP2is produced by E. coli. The target gene encodingS2-V132 is expressed with a His tag at the C-terminus. |
| Storage: | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 5375 |
| UniProt ID: | P02689 |
| Synonyms: | Myelin P2 protein; PMP2; Peripheral myelin protein 2 |



[View online »](#)