

Product datasheet for **TP726225**

PARVB Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Beta-parvin/PARVB
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-E364
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Beta-parvin/PARVBis produced by E. coli. The target gene encodingS2-E364 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	29780
UniProt ID:	Q9HBI1
Synonyms:	Beta-parvin; PARVB; CGI-56; Affixin
Protein Pathways:	Focal adhesion



[View online »](#)