

## Product datasheet for **TP726226**

### ARPC1A Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human ARPC1A Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	S2-A220
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human ARPC1A Protein is produced by E. coli. The target gene encoding S2-A220 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10552
UniProt ID:	<a href="#">Q92747</a>
Synonyms:	Actin-related protein 2/3 complex subunit 1A; ARPC1A; SOP2-like protein
Protein Pathways:	Fc gamma R-mediated phagocytosis, Pathogenic Escherichia coli infection, Regulation of actin cytoskeleton



[View online »](#)