

## Product datasheet for **TP726229**

### ZNF365 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Protein ZNF365
Species:	Human
Expression cDNA Clone or AA Sequence:	Q2-I229
Tag:	C-His
Buffer:	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	Recombinant Human Protein ZNF365is produced by E. coli. The target gene encodingQ2-I229 is expressed with a His tag at the C-terminus.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	22891
UniProt ID:	<a href="#">Q70YC4</a>
Synonyms:	Protein ZNF365; ZNF365; DISC1-binding zinc-finger protein; Protein su48



[View online »](#)